



## Validation of Quantitative Methods

Analyte	Ethanol
Unit of Measure	g/100 mL
Lowest Calibrator Concentration	0.010 g/100mL
Highest Calibrator Concentration	0.500 g/100 mL
Matrix	Blood
Analyst Performing Validation Studies	K. Baylor and M. Henry
Laboratory Director	Derek Sanders
Start Date	August 18, 2023
Completion Date	October 4, 2023
Instrument	BAC-1
Method	BAC.M

### Validation Approval

Analyst: Kayla M. Baylor 10/13/2023  
Date

Analyst: Melina Henry 10/13/2023  
Date

Laboratory Director: Derek Sanders 17 Oct 2023  
Date

Quality Manager: Ami White 20 Nov. 2023  
Date

This method validation study, even though it may depict an approval date by the crime laboratory, will not be used in casework until ANAB grants a re-accreditation to include this test method and instrument technology.

## METHOD VALIDATION PROTOCOL AND RESULTS

Analyte: Ethanol  
 Units: g/100 mL  
 Method: BAC-M  
 Instrument: BAC-1  
 SOP Reference: TOX-03-00 Validation of Quantitative Methods

Analyst: K. Baylor and M. Henry  
 Study Dates: 8/18/2023 to 10/4/2023  
 Matrix: Blood

VALIDATION EXPERIMENT		SOP CRITERIA	RESULTS	COMMENTS
1	Weight Verification	The least complex weighting scheme that minimizes 2 %RE	Unweighted: = 116.17 1/x Weighting: = 75.02 1/x <sup>2</sup> Weighting: = 69.01	1/x weighting will be used to be consistent with BAC-2 instrument. The individual calibration curves on all six days met the acceptance criteria. Also, other analytical data such as quantification values of the controls and calibrators met acceptance criteria.
1	Validation Calibration	R <sup>2</sup> value ≥ 0.99 %RE Calibrators: < ± 5% if target concentration is > 0.05 g/100mL < ± 10% if target concentration is < 0.05 g/100mL	Max %RE = 6.46	N/A
1	Linearity	95% CI of slope includes 1 95% CI of intercept includes 0	95% CI of slope = 1.0013 - 1.0075 95% CI of Intercept = -0.0016 - 0.0000	Upon inspection of the standardized residual plot, the Day 1 and Day 5, 0.500 g/100 mL calibrators were determined to be statistically significant outliers. If these data points are excluded, 95% CI of the slope contains 1. Linearity is deemed acceptable.
2	Bias & Precision	%Bias: ≤5% if target concentration is > 0.05 g/100mL ≤10% if target concentration is ≤ 0.05 g/100mL Within-Run %CV ≤10% Between-Run %CV ≤10%	Max Bias = 7.52% Max Within-Run Precision = 6.59% Max Between-Run Precision = 5.13%	Numbers in red have concentrations outside the acceptable range.
3	LOD	Signal to Noise ≥3.3 Acceptable detection and identification criteria	LOD = 0.010 g/100mL S/N = 7.9	The lowest non-zero calibrator was defined as the LOD.
4	LOQ	Bias: ≤10% Within-Run Precision: CV ≤10% Between-Run Precision: CV ≤10%	Bias = 4.45% Within-Run Precision = 3.38% Between-Run Precision = 2.25%	N/A
5	Carryover	No analyte carryover is observed in the matrix blank samples; response in blank samples is <LOQ of the method.	No significant carryover observed following samples containing ethanol at up to 1.0 g/100mL	N/A
6	Exogenous	Concentrations of analytes of interest within ± 5% of the average concentration obtained in the Bias and Precision studies	No significant interference observed.	N/A
7	Dilution	%Bias: ≤5% if diluted result is > 0.05 g/100mL ≤10% if diluted result is ≤ 0.05 g/100mL Within-Run %CV ≤10% Between-Run %CV ≤10%	It is suitable to dilute blood samples 2x. It is suitable to dilute liquid samples 20x, 50x, or 100x.	For liquid samples, numbers in red have a bias outside of 10%.
8	Autosampler Stability	Average signal (peak area, or ratio of peak area analyzer(s) compared to time 0 is within 5%	Unpunctured samples on the autosampler may be analyzed for ethanol up to 72 hours after preparation.	Two of the punctured samples did not meet the acceptance criteria when the average signal was compared to time 0. Per the SOP and this experiment, unknown samples will not be re-injected.
9	Blind Samples	N/A	All samples were within 10% of the target concentration.	Samples ran in batch BAC_20230926_MH. Blind Samples 2,3,7 and 10 were not within 10% of the target. Another licensed analyst analyzed blind samples 2,3,7 and 10 in batch BAC_20231004_KB. The samples from both batches were within 10% of each other. Therefore, the blind samples were deemed acceptable.

**Validation Study 1**

Analyte: Ethanol  
 Units: g/100 mL  
 Instrument: BAC-1

Analyte: K. Baylor and M. Henry  
 Study Dates: 8/18/2023 to 10/4/2023  
 Matrix: Blood

**STANDARD CURVE WEIGHT VERIFICATION**

Unweighted: 116.17  
 1x Weighting: 75.02  
 1x<sup>2</sup> Weighting: 69.01

Batch Name  
 BAC\_20230818\_MH

Batch Name	Cnom	Y	w	wxy	wx	wy	wx <sup>2</sup>	wy <sup>2</sup>	Found	%RE	%RE
BAC_20230818_MH		0.01	0.00518	0.01	0.05180803	0.0001	0.002684	0.01118	11.17628	11.17628	11.17628
		0.025	0.00336	0.025	0.13438403	0.000625	0.018059	0.025563	2.252028	2.252028	2.252028
		0.05	0.013793	0.05	0.27586216	0.0025	0.0761	0.050312	0.624784	0.624784	0.624784
		0.1	0.056052	0.1	0.56051693	0.01	0.314179	0.100108	0.108288	0.108288	0.108288
		0.2	0.220534	0.2	1.10267011	0.04	1.215881	0.19495	-2.52525	2.325249	
		0.4	0.901688	0.4	2.25422072	0.16	5.081511	0.396395	-0.90119	0.901186	
		0.5	2.852456	0.5	2.85245604	0.25	8.136505	0.501047	0.209432	0.209432	
BAC_20230823_KB		0.01	0.00527	0.01	0.05274451	0.0001	0.002782	0.011281	12.81452	12.81452	
		0.025	0.003343	0.025	0.13371525	0.000625	0.017788	0.025446	1.784056	1.784056	
		0.05	0.013682	0.05	0.27363668	0.0025	0.074877	0.049923	0.15384	0.15384	
		0.1	0.055092	0.1	0.55091737	0.01	0.30351	0.098429	-1.571	1.571005	
		0.2	0.223383	0.2	1.1691266	0.04	1.247494	-1.27949	1.279494		
		0.4	2.2493628	0.4	2.2493628	0.16	5.062213	0.395646	1.08856	1.08856	
		0.5	2.846594	0.5	2.84659352	0.25	8.105372	0.500092	0.18316	0.18316	
BAC_20230828_MH		0.01	0.00388	0.01	0.0388194	0.0001	0.001507	0.008845	-11.5453	11.5453	
		0.025	0.003402	0.025	0.13607572	0.000625	0.018657	0.025859	3.435765	3.435765	
		0.05	0.014065	0.05	0.28130547	0.0025	0.079133	0.051285	2.59228	2.59228	
		0.1	0.056689	0.1	0.55668769	0.01	0.309801	0.099438	-0.56158	0.561577	
		0.2	0.221827	0.2	1.10913441	0.04	1.230179	0.19608	-1.95984	1.959835	
		0.4	2.251443	0.4	2.25144287	0.16	5.068995	0.395909	-1.02267	1.022672	
		0.5	2.840772	0.5	2.84077247	0.25	8.069988	0.499003	-0.19934	0.199339	
BAC_20230918_MH		0.01	0.00053	0.01	0.05303768	0.0001	0.002813	0.011333	13.32737	13.32737	
		0.025	0.003418	0.025	0.1367267	0.000625	0.018694	0.025973	3.891282	3.891282	
		0.05	0.013717	0.05	0.27433506	0.0025	0.07526	0.050045	0.090498	0.090498	
		0.1	0.055909	0.1	0.55909245	0.01	0.312584	0.099859	-0.1409	1.40903	
		0.2	0.226404	0.2	1.13201925	0.04	1.281468	0.200084	0.041836	0.041836	
		0.4	2.262037	0.4	2.26203736	0.16	5.116813	0.397763	-0.59934	0.599337	
		0.5	2.846829	0.5	2.84682893	0.25	8.115826	0.500413	0.082531	0.082531	
BAC_20230921_MH		0.01	0.00534	0.01	0.05337596	0.0001	0.002849	0.011392	13.91914	13.91914	
		0.025	0.003399	0.025	0.13596898	0.000625	0.018488	0.02584	3.361076	3.361076	
		0.05	0.013868	0.05	0.27735199	0.0025	0.076924	0.050573	1.46028	1.46028	
		0.1	0.055872	0.1	0.55871993	0.01	0.312168	0.099794	-0.20607	0.206069	
		0.2	0.224958	0.2	1.12479023	0.04	1.265153	0.198819	-0.59046	0.590465	
		0.4	2.274661	0.4	2.27466144	0.16	5.174085	0.399971	-0.00724	0.00724	
		0.5	2.867970	0.5	2.8679704	0.25	8.340373	0.50726	1.451968	1.451968	
BAC_20230926_MH		0.01	0.000533	0.01	0.05325112	0.0001	0.002836	0.01137	13.70075	13.70075	
		0.025	0.003392	0.025	0.13569743	0.000625	0.018414	0.025793	3.171063	3.171063	
		0.05	0.013799	0.05	0.27598893	0.0025	0.07617	0.050335	0.669488	0.669488	
		0.1	0.055983	0.1	0.5598303	0.01	0.313413	0.099899	-0.01135	0.011349	
		0.2	0.226276	0.2	1.1317652	0.04	1.280017	0.199972	-0.01421	0.014206	
		0.4	2.294028	0.4	2.29402805	0.16	5.262565	0.403359	0.839731	0.839731	
		0.5	2.860276	0.5	2.8602757	0.25	8.295988	0.505914	1.182755	1.182755	
Sum		43.580382	42	15.79741	43.5803824	2.77935	89.80017				116.1658
Slope		5.71643028									

Comments: It is acceptable to use 1/x weighting because the individual calibration curves on all six days met the acceptance criteria. Also, other analytical data such as quantification values of the controls and calibrators met acceptance criteria.

Acceptance Criteria: The least complex weighting scheme that minimizes Σ|%RE|

Y	w (1/x)	wx	wy	wx^2	wy^2	Cfound	%RE	%RE
0.01	100	1	5.1808028	0.01	0.2684072	0.010367	3.668773	3.668773
0.025	40	1	5.37536113	0.025	0.7223653	0.024875	-0.49943	4.99433
0.05	20	1	5.51724323	0.05	1.521999	0.049732	-0.53546	5.35462
0.1	10	1	5.60516927	0.1	3.141792	0.099745	-0.25008	2.50056
0.2	5	1	5.6139057	0.2	6.079407	0.194999	-0.50048	2.50048
0.4	2.5	2	5.6355181	0.4	12.70378	0.397322	0.68956	6.89562
0.01	100	1	5.70491208	0.01	6.254221	0.62429	0.48584	4.8584
0.025	40	1	5.7245119	0.025	0.052745	0.010531	5.314137	5.314137
0.05	20	1	5.74460983	0.05	0.133735	0.027158	-0.39944	3.99444
0.1	10	1	5.7737356	0.1	0.273637	0.049341	-1.31748	1.31748
0.2	5	1	5.79917369	0.2	0.550917	0.090588	-1.94166	1.941659
0.4	2.5	2	5.8466328	0.4	6.237469	0.197501	-1.2493	2.49305
0.01	100	1	5.8466328	0.01	1.169133	0.396569	-0.85775	8.57752
0.025	40	1	5.8466328	0.025	2.846994	0.549694	0.2939	2.939
0.05	20	1	5.8466328	0.05	7.046664	0.25172	0.689457	6.89457
0.1	10	1	5.8466328	0.1	15.2047	0.00689	-0.91517	9.15172
0.2	5	1	5.8466328	0.2	32.4149	0.025172	0.689457	6.89457
0.4	2.5	2	5.8466328	0.4	64.8298	0.05172	0.689457	6.89457
0.01	100	1	5.8466328	0.01	0.136076	0.05172	0.689457	6.89457
0.025	40	1	5.8466328	0.025	0.281305	0.025172	0.689457	6.89457
0.05	20	1	5.8466328	0.05	0.615926	0.0136076	0.689457	6.89457
0.1	10	1	5.8466328	0.1	1.36076	0.00689	-0.91517	9.15172
0.2	5	1	5.8466328	0.2	2.81305	0.025172	0.689457	6.89457
0.4	2.5	2	5.8466328	0.4	5.6261	0.05172	0.689457	6.89457
0.01	100	1	5.8466328	0.01	0.281305	0.025172	0.689457	6.89457
0.025	40	1	5.8466328	0.025	0.615926	0.0136076	0.689457	6.89457
0.05	20	1	5.8466328	0.05	1.36076	0.00689	-0.91517	9.15172
0.1	10	1	5.8466328	0.1	2.81305	0.025172	0.689457	6.89457
0.2	5	1	5.8466328	0.2	5.6261	0.05172	0.689457	6.89457
0.4	2.5	2	5.8466328	0.4	11.2522	0.10344	0.37171	3.7171
0.01	100	1	5.8466328	0.01	0.615926	0.0136076	0.689457	6.89457
0.025	40	1	5.8466328	0.025	1.36076	0.00689	-0.91517	9.15172
0.05	20	1	5.8466328	0.05	2.81305	0.025172	0.689457	6.89457
0.1	10	1	5.8466328	0.1	5.6261	0.05172	0.689457	6.89457
0.2	5	1	5.8466328	0.2	11.2522	0.10344	0.37171	3.7171
0.4	2.5	2	5.8466328	0.4	22.5044	0.20688	0.74342	7.4342
0.01	100	1	5.8466328	0.01	0.281305	0.025172	0.689457	6.89457
0.025	40	1	5.8466328	0.025	0.615926	0.0136076	0.689457	6.89457
0.05	20	1	5.8466328	0.05	1.36076	0.00689	-0.91517	9.15172
0.1	10	1	5.8466328	0.1	2.81305	0.025172	0.689457	6.89457
0.2	5	1	5.8466328	0.2	5.6261	0.05172	0.689457	6.89457
0.4	2.5	2	5.8466328	0.4	11.2522	0.10344	0.37171	3.7171
0.01	100	1	5.8466328	0.01	0.281305	0.025172	0.689457	6.89457
0.025	40	1	5.8466328	0.025	0.615926	0.0136076	0.689457	6.89457
0.05	20	1	5.8466328	0.05	1.36076	0.00689	-0.91517	9.15172
0.1	10	1	5.8466328	0.1	2.81305	0.025172	0.689457	6.89457
0.2	5	1	5.8466328	0.2	5.6261	0.05172	0.689457	6.89457
0.4	2.5	2	5.8466328	0.4	11.2522	0.10344	0.37171	3.7171
0.01	100	1	5.8466328	0.01	0.281305	0.025172	0.689457	6.89457
0.025	40	1	5.8466328	0.025	0.615926	0.0136076	0.689457	6.89457
0.05	20	1	5.8466328	0.05	1.36076	0.00689	-0.91517	9.15172
0.1	10	1	5.8466328	0.1	2.81305	0.025172	0.689457	6.89457
0.2	5	1	5.8466328	0.2	5.6261	0.05172	0.689457	6.89457
0.4	2.5	2	5.8466328	0.4	11.2522	0.10344	0.37171	3.7171
0.01	100	1	5.8466328	0.01	0.281305	0.025172	0.689457	6.89457
0.025	40	1	5.8466328	0.025	0.615926	0.0136076	0.689457	6.89457
0.05	20	1	5.8466328	0.05	1.36076	0.00689	-0.91517	9.15172
0.1	10	1	5.8466328	0.1	2.81305	0.025172	0.689457	6.89457
0.2	5	1	5.8466328	0.2	5.6261	0.05172	0.689457	6.89457
0.4	2.5	2	5.8466328	0.4	11.2522	0.10344	0.37171	3.7171

Sum 5.69168255  
Slope -0.00718664  
Intercept 0.99898670  
R2 43.58038 1077

43.58038 72811.5 1077  
Slope 5.66668746  
Intercept -0.00622307  
R2 0.998984833

42 1277.299  
Slope 5.66668746  
Intercept -0.00622307  
R2 0.998984833


# Validation Study 1

Analyte: Ethanol  
 Units: g/100 mL  
 Instrument: BAC-1

# VALIDATION CURVE CALIBRATION

Analyst: K. Baylor and M. Henry  
 Study Dates: 8/18/2023 to 10/4/2023  
 Matrix: Blood

Batch Name	Calibration Level	Target	Calculated Result	% RE	y (Response Ratio)
BAC_20230818_MH	Level 1	0.010	0.0102	2.10	0.052
	Level 2	0.025	0.0248	0.66	0.134
	Level 3	0.050	0.0498	0.49	0.276
	Level 4	0.100	0.1000	0.01	0.561
	Level 5	0.200	0.1957	2.13	1.103
	Level 6	0.400	0.3994	0.16	2.254
	Level 7	0.500	0.5052	1.04	2.852
BAC_20230823_KB	Level 1	0.010	0.0104	3.78	0.053
	Level 2	0.025	0.0247	1.28	0.134
	Level 3	0.050	0.0494	1.14	0.274
	Level 4	0.100	0.0985	1.52	0.551
	Level 5	0.200	0.1986	0.71	1.117
	Level 6	0.400	0.3990	0.26	2.250
	Level 7	0.500	0.5046	0.91	2.847
BAC_20230828_MH	Level 1	0.010	0.0094	6.46	0.039
	Level 2	0.025	0.0265	5.90	0.136
	Level 3	0.050	0.0520	3.97	0.281
	Level 4	0.100	0.1004	0.37	0.557
	Level 5	0.200	0.1974	1.29	1.109
	Level 6	0.400	0.3981	0.48	2.251
	Level 7	0.500	0.5016	0.32	2.841
BAC_20230918_MH	Level 1	0.010	0.0102	1.88	0.053
	Level 2	0.025	0.0250	0.17	0.137
	Level 3	0.050	0.0492	1.65	0.274
	Level 4	0.100	0.0993	0.71	0.559
	Level 5	0.200	0.2001	0.06	1.132
	Level 6	0.400	0.3990	0.25	2.262
	Level 7	0.500	0.5023	0.45	2.849
BAC_20230921_MH	Level 1	0.010	0.0104	3.57	0.053
	Level 2	0.025	0.0248	0.83	0.136
	Level 3	0.050	0.0495	1.01	0.277
	Level 4	0.100	0.0987	1.32	0.559
	Level 5	0.200	0.1976	1.21	1.125
	Level 6	0.400	0.3985	0.36	2.275
	Level 7	0.500	0.5057	1.14	2.888
BAC_20230926_MH	Level 1	0.010	0.0104	3.83	0.053
	Level 2	0.025	0.0248	0.98	0.136
	Level 3	0.050	0.0492	1.60	0.276
	Level 4	0.100	0.0987	1.33	0.560
	Level 5	0.200	0.1982	0.88	1.131
	Level 6	0.400	0.4008	0.21	2.294
	Level 7	0.500	0.5030	0.60	2.880



Max %RE LOO = 6.46  
 Max %RE = 6.46

Comments: N/A

Acceptance Criteria:  
 R<sup>2</sup> value ≥ 0.99  
 %RE Calibrators:  
 < ± 5% if target concentration is > 0.05 g/100mL  
 < ± 10% if target concentration is < 0.05 g/100mL

**Validation Study 1**

Analyte: Ethanol  
 Units: g/100 mL  
 Instrument: BAC-1

**LINEARITY**

Analyst: K. Baylor and M. Henry  
 Study Dates: 8/18/2023 to 10/4/2023  
 Matrix: Blood



Batch Name	Target (x)	Calculated (y)	Predicted	Residual
BAC_20230818_MH	0.010	0.0102	-2.71	2.72
	0.025	0.0248	-2.70	2.72
	0.050	0.0498	-2.67	2.72
	0.100	0.1000	-2.62	2.72
	0.200	0.1957	-2.51	2.71
	0.400	0.3994	-2.31	2.70
	0.500	0.5052	-2.20	2.71
	0.010	0.0104	-2.71	2.72
	0.025	0.0247	-2.70	2.72
	0.050	0.0494	-2.67	2.72
BAC_20230823_KB	0.100	0.0985	-2.62	2.72
	0.200	0.1986	-2.51	2.71
	0.400	0.3990	-2.31	2.70
	0.500	0.5046	-2.20	2.71
	0.010	0.0094	-2.71	2.72
	0.025	0.0265	-2.70	2.72
	0.050	0.0520	-2.67	2.72
	0.100	0.1004	-2.62	2.72
	0.200	0.1974	-2.51	2.71
	0.400	0.3981	-2.31	2.70
BAC_20230828_MH	0.500	0.5016	-2.20	2.70
	0.010	0.0102	-2.71	2.72
	0.025	0.0250	-2.70	2.72
	0.050	0.0492	-2.67	2.72
	0.100	0.0993	-2.62	2.72
	0.200	0.2001	-2.51	2.71
	0.400	0.3990	-2.31	2.70
	0.500	0.5023	-2.20	2.70
	0.010	0.0104	-2.71	2.72
	0.025	0.0248	-2.70	2.72
BAC_20230921_MH	0.050	0.0495	-2.67	2.72
	0.100	0.0987	-2.62	2.72
	0.200	0.1976	-2.51	2.71
	0.400	0.3985	-2.31	2.70
	0.500	0.5057	-2.20	2.71
	0.010	0.0104	-2.71	2.72
	0.025	0.0248	-2.70	2.72
	0.050	0.0492	-2.67	2.72
	0.100	0.0987	-2.62	2.72
	0.200	0.1982	-2.51	2.71
BAC_20230926_MH	0.400	0.4008	-2.31	2.71
	0.500	0.5030	-2.20	2.70

Slope	1.0044
Std err in slope, $S_b$	0.0015
Degrees freedom	40
Confidence level	95%
Student t	2.0211
Confidence interval	0.003
Slope	1.004 ± 0.003
Range	1.0013 - 1.0075

NO

Intercept	-0.0008
Std err in Intercept	0.0004
Degrees freedom	40
Confidence Level	95%
Student t	2.0211
Confidence interval	0.001
Intercept	-0.001 ± 0.001
Lower	-0.0016 - 0.0000

YES



Comments: Upon inspection of the standardized residual plot, the Day 1 and Day 5, 0.500 g/100 mL calibrators were determined to be statistically significant outliers. If these data points are excluded, 95% CI of the slope contains 1. Linearity is deemed acceptable.

Acceptance Criteria:

95% CI of slope Includes 1  
 95% CI of Intercept Includes 0

**Validation Study 2**

**BIAS AND PRECISION**

Analyte: Ethanol  
 Units: g/100 mL  
 Instrument: BAC-1

Analyst: K. Baylor and M. Henry  
 Study Dates: 8/18/2023 to 10/4/2023  
 Matrix: Blood

Batch Name	Run Order	xtra low QC	Low QC	Blood QC A	Blood QC B	High QC
<i>Target Concentration (g/100 mL):</i>		0.018	0.080	0.0778	0.1976	0.400
BAC_20230818_MH	1-1	0.0188	0.0791	0.0772	0.1938	0.3995
	1-2	0.0198	0.0790	0.0770	0.2008	0.4019
	1-3	0.0189	0.0781	0.0785	0.1964	0.3989
	Mean	0.0191	0.0787	0.0776	0.1970	0.4001
	SD	0.0006	0.0005	0.0008	0.0035	0.0016
	%CV	2.91%	0.69%	1.05%	1.80%	0.41%
BAC_20230823_KB	2-1	0.0186	0.0785	0.0764	0.1966	0.3984
	2-2	0.0188	0.0785	0.0793	0.1999	0.4029
	2-3	0.0187	0.0811	0.0784	0.1986	0.4130
	Mean	0.0187	0.0794	0.0780	0.1984	0.4048
	SD	0.0001	0.0015	0.0015	0.0017	0.0075
	%CV	0.39%	1.93%	1.88%	0.85%	1.85%
BAC_20230828_MH	3-1	0.0200	0.0799	0.0779	0.1945	0.4008
	3-2	0.0205	0.0868	0.0793	0.2045	0.4246
	3-3	0.0226	0.0848	0.0793	0.2063	0.4083
	Mean	0.0210	0.0838	0.0789	0.2018	0.4112
	SD	0.0014	0.0036	0.0008	0.0064	0.0122
	%CV	6.59%	4.25%	1.02%	3.15%	2.96%
BAC_20230918_MH	4-1	0.0183	0.0794	0.0799	0.1985	0.3988
	4-2	0.0189	0.0796	0.0805	0.2094	0.4104
	4-3	0.0190	0.0817	0.0783	0.2093	0.4109
	4-4	0.0188	0.0801	--	--	0.4124
	4-5	0.0187	--	--	--	--
	Mean	0.0187	0.0802	0.0796	0.2057	0.4081
BAC_20230921_MH	5-1	0.0189	0.0782	0.0746	0.1949	0.3998
	5-2	0.0191	0.0794	0.0795	0.1948	0.4073
	5-3	0.0193	0.0806	0.0778	0.1985	0.4074
	5-4	0.0192	--	--	--	0.4075
	Mean	0.0192	0.0794	0.0773	0.1961	0.4055
	SD	0.0002	0.0012	0.0025	0.0021	0.0038
BAC_20230926_MH	6-1	--	0.0782	0.0766	0.1876	0.4002
	6-2	--	0.0808	0.0781	0.1876	0.4050
	Mean	--	0.0795	0.0774	0.1876	0.4026
	SD	--	0.0018	0.0011	0.0000	0.0034
	%CV	--	2.28%	1.37%	0.00%	0.84%
	% Bias	--	-0.60%	-0.58%	-5.06%	0.65%

Mean		0.0193	0.0802	0.0782	0.1984	0.4057
SD		0.0010	0.0023	0.0015	0.0064	0.0067
Precision (%CV)	Max Within-Run	6.59%	4.25%	3.19%	3.15%	2.96%
	Between-Run	5.13%	2.88%	1.89%	3.21%	1.66%
% Bias		7.52%	0.23%	0.40%	0.09%	1.35%

Comments: Numbers in red have concentrations outside the acceptable range.

Acceptance Criteria: %Bias: ≤5% if target concentration is > 0.05 g/100mL  
 ≤10% if target concentration is ≤ 0.05 g/100mL  
 Within-Run %CV ≤10%  
 Between-Run %CV ≤10%

### Validation Study 3

Analyte: Ethanol  
Units: g/100 mL  
Instrument: BAC-1

### Sensitivity (LOD)

Analyst: K. Baylor and M. Henry  
Study Dates: 8/18/2023 to 10/4/2023  
Matrix: Blood

Batch Name	FID-1	FID-2
Concentration (g/100 mL):	0.01	0.01
Source 1-1	8.4	81.3
Source 1-2	7.3	59.7
Source 1-3	7.2	76.5
Source 2-1	6.2	79.1
Source 2-2	9.4	76.1
Source 2-3	8.7	78.2
Source 3-1	6.5	64.5
Source 3-2	7.2	60.7
Source 3-3	9.9	60.4
Mean Signal to Noise:	7.867	70.722

See data for sample chromatogram of LOD specimen

Established LOD: 0.010 g/100mL  
Smallest Average S:N of LOD: 7.9

Comments: N/A

Acceptance Criteria:

Signal to Noise ≥ 3.3  
Acceptable detection and identification criteria



**Validation Study 4****SENSITIVITY (LOQ)**

Analyte: Ethanol  
 Units: g/100 mL  
 Instrument: BAC-1

Analyst: K. Baylor and M. Henry  
 Study Dates: 8/18/2023 to 10/4/2023  
 Matrix: Blood

Batch Name	Run Order	LOQ
<i>Target Concentration (g/100 mL):</i>		0.01
BAC_20230818_MH  <i>Within Run</i>	1-1	0.0102
	1-2	0.0104
	1-3	0.0106
	<b>Mean</b>	<b>0.0104</b>
	<b>SD</b>	<b>0.0002</b>
	<b>%CV</b>	<b>1.69%</b>
	<b>% Bias</b>	<b>3.92%</b>
BAC_20230823_KB  <i>Within Run</i>	2-1	0.0104
	2-2	0.0105
	2-3	0.0102
	<b>Mean</b>	<b>0.0104</b>
	<b>SD</b>	<b>0.0002</b>
	<b>%CV</b>	<b>1.61%</b>
	<b>% Bias</b>	<b>3.75%</b>
BAC_20230918_MH  <i>Within Run</i>	3-1	0.0102
	3-2	0.0106
	3-3	0.0109
	<b>Mean</b>	<b>0.0106</b>
	<b>SD</b>	<b>0.0004</b>
	<b>%CV</b>	<b>3.38%</b>
	<b>% Bias</b>	<b>5.69%</b>
<b>Mean</b>		<b>0.0104</b>
<b>SD</b>		<b>0.0002</b>
<b>Precision (%CV)</b>	<i>Max Within-Run</i>	<b>3.38%</b>
	<i>Between-Run</i>	<b>2.25%</b>
<b>% Bias</b>		<b>4.45%</b>

Comments: N/A

Acceptance Criteria: Bias: ≤10%  
 Within-Run Imprecision: CV ≤10%  
 Between-Run Imprecision: CV ≤10%

# Validation Study 5

Analyte: Ethanol  
Units: g/100 mL  
Instrument: BAC-1

# CARRYOVER

Analyst: K. Baylor and M. Henry  
Study Dates: 8/18/2023 to 10/4/2023  
Matrix: Blood

Average LOQ Response: \_\_\_\_\_ 7.87 \_\_\_\_\_

Sample order	Target Concentration (g/100 mL)	Results		
		Day 1	Day 2	Day 3
Date		8/18/2023	8/23/2023	8/28/23
Concentrated Sample	1.000	763.75	776.37	778.66
Blank		0.00	0.00	0.00
% of LOD Response		0.0%	0.0%	0.0%

Maximum Response in Blank: 0.0%

Comments: N/A

Acceptance Criteria:

Response in blank samples is <LOQ of the method.

## Validation Study 6

Analyte: Ethanol  
Units: g/100 mL  
Instrument: BAC-1

## EXOGENOUS SUBSTANCE INTERFERENCE

Analyst: K. Baylor and M. Henry  
Study Dates: 8/18/2023 to 10/4/2023  
Matrix: Blood

Target Concentration: 0.0802

Control Acceptance (%): 5

Batch Name: BAC\_20230818\_MH

Sample #	Compound	Amount Added	Calculated Control Concentration	% Difference from Target
1	Methanol	3.0 g/100 mL	0.0789	-1.59%
2	Isopropanol	3.0 g/100 mL	0.0777	-3.11%
3	Acetone	3.0 g/100 mL	0.0784	-2.27%
4	Ethyl Acetate	3.0 g/100 mL	0.0805	0.38%

**Conclusions:** No compound tested in this study causes significant interference with the analyte of interest.

**Comments:** N/A

**Acceptance Criteria:** Concentrations of analytes of interest within  $\pm 5\%$  of the average concentration obtained in the Bias and Precision studies

**Validation Study 7**

**DILUTION INTEGRITY**

Analyte: Ethanol  
 Units: g/100 mL  
 Instrument: BAC-1

Analyst: K. Baylor and M. Henry  
 Study Dates: 8/18/2023 to 10/4/2023  
 Matrix: Blood

Target Concentration (g/100 mL): 0.3238

Batch Name	Dilution Factor	Diluted Result (g/100mL)	Calculated Result (g/100mL)	Mean	SD	%Bias	Absolute % Bias	Within-Run Precision (per Dilution)
BAC_20230818_MH	2	0.1559	0.3118	0.3001	0.01015	-3.71	3.71	3.38
	2	0.1468	0.2936			-9.33	9.33	
	2	0.1475	0.2949			-8.92	8.92	
BAC_20230823_KB	2	0.1602	0.3204	0.3188	0.00159	-1.05	1.05	0.50
	2	0.1593	0.3186			-1.59	1.59	
	2	0.1586	0.3172			-2.03	2.03	
BAC_20230828_MH	2	0.1567	0.3133	0.3121	0.00252	-3.24	3.24	0.81
	2	0.1546	0.3092			-4.50	4.50	
	2	0.1569	0.3138			-3.08	3.08	
BAC_20230918_MH	2	0.1469	0.2937	0.3044	0.01017	-9.29	9.29	3.34
	2	0.1570	0.3140			-3.04	3.04	
	2	0.1527	0.3055			-5.65	5.65	
BAC_20230921_MH	2	0.1495	0.2990	0.3056	0.00577	-7.66	7.66	1.89
	2	0.1543	0.3087			-4.68	4.68	
	2	0.1546	0.3093			-4.49	4.49	

Dilution Factor	%Bias	Between-Run Precision (%CV)
2	-4.82	2.93

**Results:** It is suitable to dilute whole blood samples 2x prior to analysis.  
**Comments:** N/A

Overall %Bias: ≤5% if diluted result is > 0.05 g/100mL  
 Within-Run %CV ≤10%  
 Between-Run %CV ≤10%

Acceptance Criteria:

**Validation Study 7**

Analyte: Ethanol  
 Units: g/100 mL  
 Instrument: BAC-1

**DILUTION INTEGRITY**

Analyst: K. Baylor and M. Henry  
 Study Dates: 8/18/2023 to 10/4/2023  
 Matrix: Blood

Target Concentration (g/100 mL): 1.6500

Batch Name	Dilution Factor	Diluted Result (g/100mL)	Calculated Result (g/100mL)	Mean	SD	%Bias	Absolute % Bias	Within-Run Precision (per Dilution)
BAC_20230818_MH	20	0.0834	1.6690	1.6555	0.02411	1.15	1.15	1.46
	20	0.0814	1.6277			-1.35	1.35	
	20	0.0835	1.6699			1.21	1.21	
	50	0.0261	1.3049			-20.91	20.91	
	50	0.0263	1.3137	1.3086	0.00453	-20.38	20.38	0.35
	50	0.0261	1.3073			-20.77	20.77	
	100	0.0162	1.6163			-2.04	2.04	
	100	0.0160	1.6047	1.6081	0.00717	-2.75	2.75	0.45
	100	0.0160	1.6032			-2.84	2.84	
	20	0.0783	1.5668	1.5764	0.02093	-5.05	5.05	1.33
20	0.0800	1.6004			-5.34	5.34		
50	0.0319	1.5971			-3.01	3.01		
50	0.0315	1.5764			-3.21	3.21		
50	0.0317	1.5834	1.5856	0.01051	-4.46	4.46	0.66	
100	0.0166	1.6551			-4.03	4.03		
100	0.0165	1.6494	1.6547	0.00513	0.31	0.31	0.31	
100	0.0166	1.6597			-0.03	0.03		
20	0.0816	1.6329	1.6237	0.00972	0.59	0.59	0.60	
20	0.0807	1.6135			-1.03	1.03		
20	0.0812	1.6246			-2.21	2.21		
50	0.0334	1.6699			-1.54	1.54		
50	0.0336	1.6808	1.6779	0.00696	1.21	1.21	0.41	
50	0.0337	1.6828			1.87	1.87		
100	0.0177	1.7684	1.7794	0.01855	1.99	1.99	1.04	
100	0.0180	1.8009			9.14	9.14		
100	0.0177	1.7690			7.21	7.21		
20	0.0782	1.5635			-5.24	5.24		
20	0.0782	1.5632	1.5646	0.00215	-5.26	5.26	0.14	
20	0.0784	1.5671			-5.03	5.03		
50	0.0311	1.5547			-5.78	5.78		
50	0.0314	1.5691	1.5644	0.00840	-4.91	4.91	0.54	
50	0.0314	1.5694			-4.89	4.89		
100	0.0166	1.6560			0.37	0.37		
100	0.0164	1.6362	1.6501	0.01211	-0.84	0.84	0.73	
100	0.0166	1.6592			0.50	0.50		
20	0.0766	1.5325			-7.12	7.12		
20	0.0767	1.5347	1.5400	0.01121	-6.99	6.99	0.73	
20	0.0776	1.5529			-5.88	5.88		
50	0.0308	1.5393			-6.71	6.71		
50	0.0313	1.5662	1.5599	0.01831	-5.08	5.08	1.17	
50	0.0315	1.5744			-4.58	4.58		
100	0.0154	1.5424	1.5467	0.00947	-6.52	6.52	0.61	
100	0.0154	1.5401			-6.66	6.66		
100	0.0156	1.5576			-5.60	5.60		

Dilution Factor	100
%Bias	-0.13
Between-Run Precision (%CV)	4.83

Dilution Factor	50
%Bias	-6.71
Between-Run Precision (%CV)	8.29

Dilution Factor	20
%Bias	-3.51
Between-Run Precision (%CV)	2.84

Results: It is suitable to dilute whole blood samples 20x, 50x, or 100x prior to analysis.  
 Comments: Numbers in red have a bias outside of 10%.

Overall %Bias: 55% if diluted result is > 0.05 g/100mL  
 510% if diluted result is ≤ 0.05 g/100mL  
 Within-Run %CV ≤10%  
 Between-Run %CV ≤10%

Acceptance Criteria:

**Validation Study 8**

Analyte: Ethanol  
 Units: g/100 mL  
 Instrument: BAC-2

**AUTOSAMPLER STABILITY**

Analyst: K. Baylor and M. Henry  
 Study Dates: 8/17/2023 to 10/4/2023  
 Matrix: Blood

## Stability of Punctured Controls

Batch Name	Sample	Time Zero			Response Ratio	% Difference from Time Zero (Response Ratio)
		Concentration (t0)	Response Ratio	Concentration (t0 72 hours later)		
BAC_20230919_KB and BAC_20230922_KB	Low QC	0.0790	0.4526	0.08390	0.4760	-5.19%
		0.0782	0.4482	0.08277	0.4695	-4.75%
		0.0796	0.4561	0.08414	0.4774	-4.68%
BAC_20230922_KB and BAC_20230925_KB	High QC	0.4048	2.3132	0.41909	2.40	-3.83%
		0.4041	2.3088	0.42150	2.42	-4.62%
		0.4017	2.2951	0.41625	2.39	-3.93%

## Stability of Unpunctured Controls

Batch Name	Sample	Time Zero			Response Ratio	% Difference from Time Zero (Response Ratio)
		Concentration (t0)	Response Ratio	Concentration (t72)		
BAC_20230919_KB and BAC_20230922_KB	Low QC	0.0790	0.4526	0.08230	0.4669	-3.16%
		0.0782	0.4482	0.08174	0.4636	-3.44%
		0.0796	0.4561	0.08228	0.4668	-2.35%
BAC_20230922_KB and BAC_20230925_KB	High QC	0.4048	2.3132	0.41729	2.3913	-3.38%
		0.4041	2.3088	0.41193	2.3605	-2.24%
		0.4017	2.2951	0.40907	2.3441	-2.14%

**Results:** Unpunctured samples on the autosampler may be analyzed for ethanol up to 72 hours after preparation.

**Comments:** One of the punctured samples did not meet the acceptance criteria when the average signal was compared to time 0. Per the SOP, unknown samples will not be re-injected.

Acceptance Criteria:

Average signal (peak area, or ratio of peak area analyte/IS) compared to time 0 is within 5%.

# Validation Study 9

Analyte: Ethanol  
Units: g/100 mL  
Instrument: BAC-1

# BLIND SAMPLES

Analyst: K. Baylor and M. Henry  
Study Dates: 8/18/2023 to 10/4/2023  
Matrix: Blood

Batch Name BAC\_20230926\_MH

Sample ID	Expected Results	Reported Results	% Difference	Pass/Fail
Blind Sample 1	0.0789	0.0835	5.79%	Pass
Blind Sample 2	0.2084	0.2185	4.82%	Pass
Blind Sample 3	0.0219	0.0232	6.03%	Pass
Blind Sample 4	Neg	Neg	-	Pass
Blind Sample 5	0.1578	0.1707	8.15%	Pass
Blind Sample 6	0.3945	0.4199	6.44%	Pass
Blind Sample 7	0.0471	0.0456	-3.18%	Pass
Blind Sample 8	0.1775	0.1929	8.65%	Pass
Blind Sample 9	Neg	Neg	-	Pass
Blind Sample 10	0.0682	0.0695	1.79%	Pass

**Comments:** Samples ran in batch BAC\_20230926\_MH. Blind Samples 2,3,7 and 10 were not within 10% of the target. Another licensed analyst analyzed blind samples 2,3,7 and 10 in batch BAC\_20231004\_KB. The samples from both batches were within 10% of each other. Therefore, the blind samples were deemed acceptable.

Acceptance Criteria:

Blind Samples are within 10% of the expected.

## SUMMARY OF VALIDATION PERFORMANCE

Analyte: Ethanol  
Units: g/100 mL  
Instrument: BAC-1

Analyst: K. Baylor and M. Henry  
Study Dates: 8/18/2023 to 10/4/2023  
Matrix: Blood

The intent of this summary is to capture and document important information about the performance of this method outside the required measurements for validation.

Failed Runs (include dates/reasons):

Date	Reason
N/A	N/A

Deviations from SOP:

N/A

Other Observations:

Reagent Verified in Validation:  
Internal Standard: R-23-5002

Sample Preparation Steps:

Refer to TOX-02-00 Ethanol Analysis Using Headspace Gas Chromatography

Location of Raw Data:

BAC1 Validation 2023 binder.

Recommended Maximum Run Length (Unknown Samples):

30

Conclusion:

This method is fit for use on casework for ethanol confirmation analysis in blood and liquid samples.