



Verification of Quantitative Methods

Purpose	To show that the instrument is fit for use after completion of annual preventative maintenance.
Analyte	Ethanol
Unit of Measure	g/100 mL
Lowest Calibrator Concentration	0.010 g/100mL
Highest Calibrator Concentration	0.500 g/100 mL
Primary Matrix	Blood
Secondary Matrices	Alcoholic Beverages and Other Liquid Specimens
Analyst Performing Verification Study	Melissa Henry
Reviewer	Kayla Baylor
Start Date	January 30, 2024
Completion Date	January 30, 2024
Instrument	BAC-1
Method	BAC.M

Validation Approval

Analyst: _____ Date

Reviewer: _____ Date

Verification Study		BIAS AND PRECISION			
Analyte: Ethanol		Analyst: Melissa Henry			
Units: g/100 mL		Study Dates: 1/30/2024			
Instrument: BAC-1	FID1	Matrix: Blood			
Batch Name	Run Order	Low QC	Blood QC A	Blood QC B	High QC
<i>Target Concentration (g/100mL):</i>		0.080	0.0778	0.1976	0.400
BAC_20240130_MH <i>Within Run</i>	1	0.0764	0.0763	0.1926	0.4033
	2	0.0768	0.0763	0.1942	0.4085
	3	0.0774	0.0765	0.1966	0.4159
	Mean	0.0769	0.0764	0.1945	0.4092
	SD	0.0005	0.0001	0.0020	0.0063
	%CV	0.65%	0.15%	1.04%	1.55%
	% Bias	-3.92%	-1.84%	-1.59%	2.31%
Comments:					
Acceptance Criteria: %Bias: ≤5% if target concentraton is > 0.05 g/100mL ≤10% if target concentraton is ≤ 0.05 g/100mL Within-Run %CV ≤10%					

Verification Study		BIAS AND PRECISION			
Analyte: Ethanol		Analyst: Melissa Henry			
Units: g/100 mL		Study Dates: 1/30/2024			
Instrument: BAC-1	FID2	Matrix: Blood			
Batch Name	Run Order	Low QC	Blood QC A	Blood QC B	High QC
<i>Target Concentration (g/100mL):</i>		0.080	0.0778	0.1976	0.400
BAC_20240130_MH <i>Within Run</i>	1	0.0757	0.0759	0.1925	0.4058
	2	0.0764	0.0759	0.1944	0.4112
	3	0.0775	0.0762	0.1971	0.4190
	Mean	0.0765	0.0760	0.1947	0.4120
	SD	0.0009	0.0002	0.0023	0.0066
	%CV	1.19%	0.23%	1.19%	1.61%
	% Bias	-4.33%	-2.31%	-1.48%	3.00%
Comments:					
Acceptance Criteria: %Bias: ≤5% if target concentraton is > 0.05 g/100mL ≤10% if target concentraton is ≤ 0.05 g/100mL Within-Run %CV ≤10%					

SUMMARY OF VALIDATION PERFORMANCE

Analyte: **Ethanol**
Units: **g/100 mL**
Instrument: **BAC-1**

Analyst: **Melissa Henry**
Study Dates: **1/30/2024**
Matrix: **Blood**

Failed Runs (include dates/reasons):

Date	Reason
N/A	N/A

Deviations from SOP: N/A

Notes: N/A

Conclusion: **BAC-1 is fit for analysis of ethanol casework.**